

# B-BOX™ Blue Light LED epi-illuminator

Catalogue No.: VE0100

## Operator's Manual

2022 ver. 1.4.0



## Product check list

- VE0100 B-BOX™ Blue Light LED epi-illuminator
- VE0101 Adapter, AC 100-240V, 50/ 60 Hz, 1.8m Power Cord
- VE0102 Multi-angle Filter Plate (0°, 30°, 60° and 90°) for B-BOX™, Acrylic,  
Amber Color, 154 x 194 (mm)
- VE0103 B-BOX™ shades

Operator's Manual and Warranty Certificate

2013 Product of SMOBiO, Inc. All rights reserved, no part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, in whole or part, without the prior written permission of SMOBiO, Inc.

Patents Pending



B-BOX™,



SMOBiO logos are registered trademarks of


SMOBiO, Inc.

## **Contents**

- A. Introduction
- B. Important Safety Information
- C. Description and Specifications
- D. Installation and Operation
- E. Warranty information
- F. Troubleshooting
- G. Information on ordering
- H. Related Products

## A. Introduction

The SMOBiO B-BOX™ Blue Light LED epi-illuminator is designed to meet the needs of a light and compact bench top illumination unit for visualizing electrophoresis. B-BOX™ features light-emitting diodes or LED light modules which emit high intensity blue light of 470nm wavelength, while emitting a minute amount of heat. B-BOX is compatible with non-Ethidium bromide (EtBr), fluorescent DNA and protein dyes. B-Box offers an unparalleled level of detection sensitivity and protection of DNA and protein samples from UV damage.




## B. Important Safety Information

- B. Important safety information
- \* Staring directly into the LED without protective Eyeware or a filter plate may cause damage to the retina; proper safety measures or protocols must be always be followed/applied.
- \* The working area is water tight, but B-BOX™ is not waterproof. Do not submerge or pour excessive amounts of liquid onto the B-BOX™ epi-illuminator.
- \* High intensity, low-heat LED modules emit only a small amount of heat during operation. However, as a precaution, the B-BOX epi-illuminator should be placed in a well-ventilated area with enough vertical clearance for airspace at the baseplate while (in use or being used).

## C. Description and Specifications


- Overall dimension (mm): 201.4 x 200 x 38 (D x W x H)
- Viewing area (mm): 158 x 96 (D x W)
- LED wavelength (nm): 470

- Number of LED Units: 72 PCS Super Flux LED
- LED life span (470 nm), up to 50,000 hours
- Power : DC 12 Volt, 0.72 Amp
- Electrical Requirements: AC 100~240V, 50/60 Hz (Adapter)
- Weight (kg): 0.95 (Net Weight)



## D. Installation and Operation

Carefully unpack the B-Box Blue Light LED epi-illuminator. The epi-illuminator should be placed on a dry horizontal surface following the important safety measures previously outlined in section B. Before operation the B-BOX Blue Light Epi-illuminator must be plugged into a compatible power source as specified in section C. When operating the epi-illuminator, place the gel on top of the working area after all excess liquids, staining, or destaining solutions have been sufficiently drained. The included **B-BOX shades** (protective eye goggles/shades) or the filter plate must be used when the LED module set is in use. Gel cutting may be performed directly on top of the glass surfaced working area. The filter plate is available as a four-angled design for storage at 0 (degrees), visualization at 30 and 90 degrees, and for gel cutting at 60 degrees)



### Caution!!!

- Do not place any heavy objects on the filter plate or bend it forcefully.
- Avoid contact with rough surfaces.
- Avoid using organic solvents such as isopropanol or ethanol while cleaning the illuminator.
- For desired cleaning results use a mild detergent and deionized water.
- Do not submerge the illuminator into any liquid solution.
- Do not exert excessive pressure on the working area or any part of the B-BOX unit.
- The B-Box Blue Light LED epi-illuminator is compatible with the following fluorescent dyes:
- FluoroalStain Fluorescent Staining dyes+ (we may say signifies a SMOBiO product or an SMOBiO product)
  
- FluoroalStain Fluorescent Staining Dyes from SMOBiO
- SYBR Green I/II
- GFP
- SYBR Safe/Gold
- SafeView
- Midori Green Advanced
- Gel Red/Green
- HealthView
- SERVA DNA Stain Clear G
- EZ-Vision

## E. Warranty information

The B-BOX™ Blue Light LED epi-illuminator is warranted against defects and workmanship for one year. Any defect which occurred within the warranty period will be repaired. Any defective part will be replaced by SMOBiO without charge.

The following circumstances the warranty will not be covered under warranty:

- \*The warranty was not registered or no proof of purchase exists.
- \*Improper user error.
- \*Unauthorized modification and replacement of parts by an unauthorized party.
- \*Damage caused by substituting alternative parts supplied by an unauthorized party.
- \*Damage caused by accident or misuse. (the same as 2nd one)
- \*Damage caused by disaster. (fire?)
- \*Corrosion caused by the use of improper solvents or samples.

## F. Trouble Shooting

G. Listed below are possible problems with recommended solutions. If you cannot resolve the problem, please call or e-mail SMOBiO, Inc. support@smobio.com or contact your Distributor.

| <b>Problem</b>              | <b>Causes</b>   | <b>Solutions</b>  |
|-----------------------------|---|---|
| Leakage of the working area | <ul style="list-style-type: none"><li>• Excessive pressure being applied to the surface or improper use of organic solvent in cleaning.</li></ul>   | <ul style="list-style-type: none"><li>• Contact your local distributor for service.</li></ul>   |
| Malfunctioning LED modules  | <ul style="list-style-type: none"><li>• LED Modules are not working correctly.</li></ul>  | <ul style="list-style-type: none"><li>• Contact your local distributor for service. Do not attempt to replace parts of the B-Box Blue Light LED epi-illuminator. Replacing any parts will void the warranty of the product.</li></ul>   |
| High background             | <ul style="list-style-type: none"><li>• Increased background signal. Contaminants within the electrophoresis running buffer may cause local deposition of staining material on the gel resulting in undesired background signals.</li></ul> | Sufficiently drain DNA or protein gel prior to visualization to reduce background signal.<br>Thoroughly drain off running buffer, or give it a quick rinse with deionized water. Refer to the manufacturer's instruction manual for the dye being used. <ul style="list-style-type: none"><li>•</li></ul> |


## H. Ordering Information

|        |  |
|--------|--|
| VE0100 | B-BOX™ Blue Light LED epi-illuminator, DC 12 Volt, 0.72 Amp<br>Complete Set Includes: VE0101, VE0102, VE0103 and Operator's Manual |
| VE0101 | Adapter, AC 100-240V, 50/ 60 Hz, 1.8m Power Cord   |
| VE0102 | Multi-angle Filter Plate (0°, 30°, 60° and 90°) for B-BOX™, Acrylic, Amber Color, 154 x 194 (mm)                                   |
| VE0103 | Goggles, Amber Color   |

## I. Related Products

|        |   |
|--------|---|
| DM1160 | FluoroBand 50 bp Fluorescent DNA Ladder, 500 µl                               |
| DM2160 | FluoroBand 100 bp Fluorescent DNA Ladder, 500 µl                              |
| DM2360 | FluoroBand 100 bp+3K Fluorescent DNA Ladder, 500 µl                           |
| DM3160 | FluoroBand 1 KB (0.25-10 kb) Fluorescent DNA Ladder, 500 µl                   |
| DM3260 | FluoroBand 1 KB Plus (0.1-10 kb) Fluorescent DNA Ladder, 500 µl               |
| DM4160 | FluoroBand XL 25 kb Fluorescent DNA Ladder, Broad Range (up to 25 kb), 500 µl |
| DL5000 | FluoroDye DNA Fluorescent Loading Dye (Green, 6X), 1 ml                       |
| DL5001 | FluoroDye DNA Fluorescent Loading Dye (Green, 6X), 1 ml x 5                   |
| DS1000 | FluoroStain DNA Fluorescent Staining Dye (Green, 10,000X), 500 µl             |
| DS1001 | FluoroStain DNA Fluorescent Staining Dye (Green, 10,000X), 500 µl x 5         |
| PS1000 | FluoroStain Protein Fluorescent Staining Dye (Red, 1,000X), 1 ml              |
| PS1001 | FluoroStain Protein Fluorescent Staining Dye (Red, 1,000X), 1 ml x 5          |
| TP1260 | ExcelTaq Fluorescent 5X PCR Master Mix, 200 Rxn                               |

## Distributed by:



The latest version of the manual can be downloaded from [www.smobio.com/shop](http://www.smobio.com/shop).

Manufactured by SMOBiO, Inc.

E-mail: [smobio@smobio.com](mailto:smobio@smobio.com)



[www.smobio.com](http://www.smobio.com)