

pSK-BNT-Swt

LOCUS pSK-BNT-Swt 7197 bp DNA circular SYN 14-FEB-2022
SOURCE synthetic DNA construct.
ORGANISM synthetic DNA construct
REFERENCE 1 (bases 1 to 7197)
AUTHORS Demo User
TITLE Direct Submission
JOURNAL Exported Friday, Jan 21, 2022 from SnapGene Viewer 5.1.7 to Vector
NTI(R) format
<https://www.snapgene.com>
COMMENT Storage Type: Basic
Original Author:
Demo User
COMMENT This file is created by Vector NTI
<http://www.invitrogen.com/>
COMMENT ORIGDB|GenBank
COMMENT VNTDATE|-14264406|
COMMENT VNTDBDATE|-14264406|
COMMENT LSOWNER|
COMMENT VNTNAME|pSK-BNT-Swt|
COMMENT VNTAUTHORNAME|Demo User|
FEATURES Location/Qualifiers
misc_feature complement(263..569)
/vntifkey="21"
/label=f1
misc_feature 754..772
/vntifkey="21"
/label=T7\promoter
misc_feature 771..823
/vntifkey="21"
/label=5'\UTR
misc_feature 773..5070
/vntifkey="21"
/label=Cassette\A
primer 794..823
/vntifkey="27"
/label=PBNT5end
/note="sequence: TCCCCACAGACTCAGAGAGAACCCGCCACC"
primer complement(794..823)
/vntifkey="27"
/label=MBNT5end
/note="sequence: GGTGGCGGGTTCTCTCTGAGTCTGTGGGGA"
CDS 824..4642
/codon_start=1
/vntifkey="4"
/label=S\protein
primer 2060..2080
/vntifkey="27"
/label=PBNT1292
/note="sequence: ggacagacaggcaagatcgcc"
primer 2132..2161
/vntifkey="27"
/label=PBNT1ig1
/note="sequence: AACAGCAACAACCTGGACTCCAAAGTCGGC"
primer complement(2132..2161)
/vntifkey="27"

pSK-BNT-Swt

```

/misc_feature /label=MBNT1ig1
                /note="sequence: GCCGACTTTGGAGTCCAGGTTGTTGCTGTT"
                2177..2179
                /vntifkey="21"
/misc_feature /label=L452
                2273..2275
                /vntifkey="21"
                /label=E484
/misc_feature 2324..2326
                /vntifkey="21"
                /label=N501
/misc_feature 2663..2664
                /vntifkey="21"
                /label=D614
/misc_feature 2864..2866
                /vntifkey="21"
                /label=P681
/primer        3461..3490
                /vntifkey="27"
                /label=PBNT1ig2
                /note="sequence: GGCACAATCACAAGCGGCTGGACATTTGGA"
/primer        complement(3461..3490)
                /vntifkey="27"
                /label=MBNT1ig2
                /note="sequence: TCCAAATGTCCAGCCGCTTGTGATTGTGCC"
/primer        3653..3691
                /vntifkey="27"
                /label=P-BNT950N
                /note="sequence: gccctgggaaagctgcagaacgtggtcaaccagaatgcc"
/misc_feature 3671..3673
                /vntifkey="21"
                /label=D950
/primer        4139..4159
                /vntifkey="27"
                /label=PBNT3387
                /note="sequence: cagcggaaacttctacgagccc"
/primer        4652..4681
                /vntifkey="27"
                /label=PBNT3end
                /note="sequence: GAGCTGGTACTGCATGCACGCAATGCTAGC"
/primer        complement(4652..4681)
                /vntifkey="27"
                /label=MBNT3end
                /note="sequence: GCTAGCATTGCGTGCATGCAGTACCAGCTC"
/misc_feature 4655..4937
                /vntifkey="21"
                /label=3'\UTR
/misc_feature 4944..5054
                /vntifkey="21"
                /label=polyA
/primer        complement(5136..5158)
                /vntifkey="27"
                /label=T3P
                /note="sequence: cgcaattaaccctcactaaaggg"
/misc_feature  complement(5136..5155)

```

pSK-BNT-Swt

```

/vnt ifkey="21"
/label=T3\P
primer complement(5460..5481)
/vnt ifkey="27"
/label=PT6265\seq
primer /note="sequence: cgtattaccgcctttgagtgag"
5514..5537
/vnt ifkey="27"
/label=P-DL-1
misc_feature /note="sequence: agaacatgtgagcaaaaggccagc"
complement(5522..6189)
/vnt ifkey="21"
/label=pUC\ori
primer complement(6193..6214)
/vnt ifkey="27"
/label=M-Pica-3561
misc_feature /note="sequence: ttcggtccactgagcgtcagac"
complement(6340..7197)
/vnt ifkey="21"
/label=Ampicillin

```

BASE COUNT 1843 a 2040 c 1799 g 1515 t
ORIGIN

```

1  acttttcctt  tttcaatatt  attgaagcat  ttatcagggt  tattgtctca  tgagcggata
61  catatttgaa  tgtatttaga  aaaataaaca  aataggggtt  cgcgcacat  ttccccgaaa
121 agtgccacct  aaattgtaag  cgttaatatt  ttgttaaaat  tcgcgttaa  tttttgttaa
181 atcagctcat  tttttaacca  ataggccgaa  atcggcaaaa  tcccttata  atcaaaagaa
241 tagaccgaga  tagggttgag  tgttgttcca  gtttggaa  agagtccact  attaaagaac
301 gtggactcca  acgtcaaagg  gcgaaaaacc  gtctatcagg  gcgatggccc  actacgtgaa
361 ccataccct aatcaagttt  ttgggggtcg  aggtgccgta  aagcaactaa  tcggaacct
421 aaagggagcc  cccgatttag  agcttgacgg  ggaaagccgg  cgaacgtggc  gagaaaggaa
481 gggaaagaa  cgaaaggagc  gggcgctagg  gcgctggcaa  gtgtagcgg  cacgctgcgc
541 gtaaccacca  caccgcgcgc  gcttaatgcg  ccgctacagg  gcgctgcca  ttccgcatc
601 aggctgcgca  actgttggga  agggcgatcg  gtgcgggct  ctccgctatt  acgccagctg
661 gcgaaagggg  gatgtgctgc  aaggcgatta  agttgggtaa  cggcagggtt  ttcccagtca
721 cgacgttgta  aaacgacggc  cagtgagcgc  gcgtaatacg  actcactata  ggaataaact
781 agtattcttc  tggteccac  agactcagag  agaaccgcgc  accatgttcg  tgttctctgt
841 gctgctgct  ctgggtgtcc  gccagtgtgt  gaacctgacc  accagaacac  agctgcctcc
901 agcctacacc  aacagcttta  ccagaggcgt  gtactacccc  gacaagggtg  tcagatccag
961 cgtgctgca  tctaccagg  acctgttct  gcctttcttc  agcaacgtga  cctggttcca
1021 cgccatccac  gtgtccggca  ccaatggcac  caagagattc  gacaaccccg  tctgctcctt
1081 caacgacggg  gtgtactttg  ccagcaccga  gaagtccaac  atcatcagag  gctggatctt
1141 cggcaccaca  ctggacagca  agaccagag  cctgctgac  gtgaacaacg  ccaccaacgt
1201 ggtcatcaaa  gtgtgagagt  tccagttctg  caacgacccc  ttcttgggcg  tctactacca
1261 caagaacaac  aagagctgga  tggaaagcga  gttccgggtg  tacagcagcg  ccaacaactg
1321 cacttcgag  tacgtgtccc  agccttctct  gatggacctg  gaaggcaagc  agggcaactt
1381 caagaacctg  cgcgagttcg  tgtttaagaa  catcgacggc  tacttcaaga  tctacagcaa
1441 gcacaccct  atcaacctcg  tgcgggatct  gcctcagggc  ttctctgctc  tggaaacctt
1501 ggtggatctg  cccatcgga  tcaacatcac  ccggtttcag  aactgtctgg  cctgacacag
1561 aagctacctg  acacctggcg  atagcagcag  cggatggaca  gctgggtgccc  ccgcttacta
1621 tgtgggctac  ctgcagccta  gaaccttct  gctgaagtac  aacgagaacg  gcaccatcac
1681 cgacgcctg  gattgtgctc  tggatcctct  gagcgagaca  aagtgcaccc  tgaagtcttt
1741 caccgtggaa  aagggcatct  accagaccag  caacttccgg  gtgcagccca  ccgaatccat
1801 cgtgcggttc  cccaataca  ccaatctgtg  ccccttccgg  gaggtgttca  atgccaccag
1861 attgectct  gtgtacgct  ggaaccgga  gcggatcagc  aattgctgtg  ccgactactc
1921 cgtgctgtac  aactccgcca  gcttcagcac  ctccaagtgc  tacggcgtgt  ccctaccaa

```

pSK-BNT-Swt

1981 gctgaacgac ctgtgcttca caaacgtgta cgccgacagc ttcgtgatcc ggggagatga
 2041 agtgcggcag attgcccctg gacagacagg caagatcgcc gactacaact acaagctgcc
 2101 cgacgacttc accggctgtg tgattgcctg gaacagcaac aacctggact ccaaagtctg
 2161 cggcaactac aattacctgt accggctgtt ccggaagtcc aatctgaagc ccttcgagcg
 2221 ggacatctcc accgagatct atcaggcccg cagcaccctt tgaacggcg tggaaaggctt
 2281 caactgctac tcccactgc agtcctacgg ctttcagccc acaaattggcg tgggctatca
 2341 gccctacaga gtggtggtgc tgagcttcca actgctgcat gccctgcca cagtgtgctg
 2401 ccctaagaaa agcaccaatc tcgtgaagaa caaatgctg aacttcaact tcaacggcct
 2461 gaccggcacc ggcgtgctga cagagagcaa caagaagtcc ctgccattcc agcagtttgg
 2521 ccgggataat gccgatacca cagacgccgt tagagatccc cagacactgg aaatcctgga
 2581 caccaccctt tgcagcttcc ggggagtgtc tgtgatcacc cctggcacca acaccagcaa
 2641 tcaggtggca gtgctgtacc aggacgtgaa ctgtaccgaa gtgcccgtgg ccattcacgc
 2701 cgatcagctg acacctacat ggcgggtgta ctccaccggc agcaatgtgt ttcagaccag
 2761 agccggctgt ctgatcggag ccgagcacgt gaacaatagc tacgagtggc acatccccat
 2821 cggcgttggg atctgcgcca gctaccagac acagacaaac agccctcgga gagccagaag
 2881 cgtggccagc cagagcatca ttgcctacac aatgtctctg ggcgcgcgaga acagcgtggc
 2941 ctactccaac aactctatcg ctatccccac caacttcacc atcagcgtga ccacagagat
 3001 cctgctgtgt tccatgacca agaccagcgt ggactgcacc atgtacatct gcggcgattc
 3061 caccgagtgc tccaacctgc tgcctgcagta cggcagcttc tgcaccagc tgaatagagc
 3121 cctgacaggg atcggcgtgg aacaggacaa gaacacccaa gaggtgttcc cccaagtga
 3181 gcagatctac aagaccctc ctatcaagga cttcggcggc tccaatttca gccagattct
 3241 gccgatcct agcaagccca gcaagcggag ctctatcgag gacctgctgt tcaacaaagt
 3301 gacactggcc gacgccggct tcatcaagca gtatggcgat tgtctgggcg acattgccgc
 3361 cagggatctg atttgcgccc agaagtttaa cggactgaca gtgctgcctc ctctgctgac
 3421 cgatgagatg atcggcccgt acacatctgc cctgctggcc ggcacaatca caagcggctg
 3481 gacatttggg gcagggcggc ctctgcagat cccctttgct atgcagatgg cctaccggtt
 3541 caacggcacc ggagtgacc agaatgtgct gtacgagaac cagaagctga tgcaccaaca
 3601 gttcaacagc gccatcgcca agatccagga cagcctgagc agcacagcaa gcgccctggg
 3661 aaagtgcag gacgtggtca accagaatgc ccaggcactg aacaccctgg tcaagcagct
 3721 gtctccaac ttcggcgcca tcagctctgt gctgaacgat atcctgagca gactggacc
 3781 tctgaggcc gaggtgcaga tcgacagact gatcacaggc agactgcaga gcctccagac
 3841 atactgacc cagcagctga tcagagccgc cgagattaga gcctctgcca atctggccgc
 3901 caccaagatg tctgagtgtg tgcctggcca gagcaagaga gtggactttt gcggcaaggg
 3961 ctaccacctg atgagcttcc ctgagctgac cctcacggc gtggtgtttc tgcagctgac
 4021 atatgtgccc gctcaagaga agaatttcac caccgctcca gccatctgcc acgacggcaa
 4081 agcccacttt cctagagaag gctgttctgt gtccaacggc acccattggt tctgacaca
 4141 gcggaacttc tacgagcccc agatcctcac caccgacaac acctctgtgt ctggcaactg
 4201 cgactctgtg atcggcattg tgaacaatac cgtgtacgac cctctgcagc ccgagctgga
 4261 cagcttcaaa gaggaactgg acaagtactt taagaaccac acaagccccg acgtggacct
 4321 gggcgataat agcggaaatc atgccagcgt cgtgaacatc cagaaagaga tgcaccggct
 4381 gaacgaggtg gccaaagaatc tgaacgagag cctgatcgac ctgcaagaac tggggaagta
 4441 cgagcagtac atcaagtggc cctggctacat ctggctgggc tttatcgccg gactgattgc
 4501 catcgtgatg gtcacaatca tgcctgtgtg catgaccagc tgcctgtagct gcctgaaggg
 4561 ctgttctgag tctggcagct gctgcaagtt cgacgaggac gattctgagc ccgtgctgaa
 4621 gggcgtgaaa ctgcaactaca catgatgact cgagctggta ctgcatgcac gcaatgctag
 4681 ctgccccctt cccgtctctgg gtacccccgag tctccccga cctcgggtcc caggtatgct
 4741 cccacctcca cctgccccac tcaccacctc tgcctagtcc agacacctcc caagcacgca
 4801 gcaatgcagc tcaaaacgct tagcctagcc acacccccac gggaaacagc agtgattaac
 4861 ctttagcaat aaacgaaagt ttaactaagc tatactaacc ccagggttgg tcaatttctg
 4921 gccagccaca ccttggagct agcaaaaaaaaa aaaaaaaaaa aaaaaaaaaa aaaaaaaaaa
 4981 actaaaaaaaa aaaaaaaaaa aaaaaaaaaa aaaaaaaaaa aaaaaaaaaa aaaaaaaaaa
 5041 aaaaaaaaaa aaaagcttga tctgaattc ctgcagcccc ggggatccac tagttctaga
 5101 gggcggcca ccgcggtgga gctccagctt ttgttccctt tagtgagggt taattgcgcg
 5161 ctggcgtaa tcatggtcat agctgtttcc tgtgtgaaat tgttatccgc tcacaattcc
 5221 acacaacata cgagccggaa gcataaagtg taaagcctgg ggtgcctaatt gagtgagcta

pSK-BNT-Swt

5281 actcacatta attgcgttgc gctcactgcc cgctttccag tcgggaaacc tgtcgtgcca
 5341 gctgcattaa tgaatcggcc aacgcgcggg gagaggcggg ttgcgtattg ggcgctcttc
 5401 cgcttctctg ctactgact cgctgcgctc ggtcgttcgg ctgcggcgag cggtatcagc
 5461 tactcaaaag gcgtaatac ggttatccac agaatacagg gataacgcag gaaagaacat
 5521 gtgagcaaaa ggccagcaaa aggccaggaa ccgtaaaaag gccgcgttgc tggcgTTTTT
 5581 ccataggctc cgccccctg acgagcatca caaaaatcga cgctcaagtc agaggTggcg
 5641 aaaccgaca ggactataaa gataccaggc gtttccccct ggaagctccc tcgtgcgctc
 5701 tctgttccg acctgcccgc ttaccggata cctgtccgcc tttctccctt cgggaagcgt
 5761 ggcgctttct catagctcac gctgtaggta tctcagttcg gtgtaggctg ttgcctccaa
 5821 gctgggctgt gtgcacgaac cccccgttca gcccgaccgc tgcgccttat ccgtaacta
 5881 tcgtcttgag tccaaccggg taagacacga cttatcgcca ctggcagcag ccaactggtaa
 5941 caggattagc agagcgaggT atgtaggcgg tgctacagag ttcttgaagt ggtggcctaa
 6001 ctacggctac actagaagga cagtatttgg tatctgcgct ctgctgaagc cagttacctt
 6061 cggaaaaaga gttggtagct cttgatccgg caaacaacc acccgtggta gcggtggttt
 6121 ttttgtttgc aagcagcaga ttacgcgcag aaaaaagga tctcaagaag atcctttgat
 6181 ctttctacg gggctctgacg ctcagtggaa cgaaaactca cgttaaggga ttttggctat
 6241 gagattatca aaaaggatct tcacctagat ctttttaaat taaaaatgaa gttttaaatc
 6301 aatctaaagt atatatgagt aaacttggtc tgacagttac caatgcttaa tcagtgaggc
 6361 acctatctca gcgatctgtc ttttctgtc atccatagtt gctgactcc ccgtcgtgta
 6421 gataactacg atacgggagg gcttaccatc tggccccagt gctgcaatga taccgcgaga
 6481 cccacgctca ccggctccag atttatcagc aataaaccag ccagccggaa gggccgagcg
 6541 cagaagtggT cctgcaactt tatccgcctc catccagtct ataatgtgtt gccgggaagc
 6601 tagagtaagt agttccccag ttaatagttt gcgcaacgtt gttgccattg ctacaggcat
 6661 cgtgggtgta cgctcgtcgt ttgggatggc ttcatcagc tccggttccc aacgatcaag
 6721 gcgagttaca tgatccccc tgttgtgcaa aaaagcgggtt agctccttcg gtccctccgat
 6781 cgttgtcaga agtaagtTgg ccgcagtgtt atcactcatg gttatggcag cactgcataa
 6841 ttctcttact gtcatgccaT ccgtaagatg cttttctgtg actgggtgagt actcaaccaa
 6901 gtcattctga gaatagtTga tgcggcgacc gagttgctct tgcocggcgt caatacggga
 6961 taatacccgC ccacatagca gaactttaa agtgcctcat attggaaaac gttcttcggg
 7021 gcgaaaactc tcaaggatct taccgctgtt gagatccagt tcgatgtaac ccaactcgtgc
 7081 acccaactga tcttcagcat cttttacttt caccagcgtt tctgggtgag caaaaacagg
 7141 aaggcaaaat gccgcaaaa agggaataag ggcgacacgg aaatgttTgaa tactcat

//